



## VENTILATOR BMC Y SERIES BPAP SYSTEM

**€3,295.00**

**SKU: N/A**

**Categories:** [Bioide](#), [Machines](#), [Ventilators](#)

### PRODUCT DESCRIPTION

Y-30T

Dimensions: 170 mm × 180 mm × 118 mm 290 mm × 180 mm × 134 mm (with the humidifier)

Weight: 1.5 kg  
2.5 kg (with the humidifier)

Water capacity: 350 mL at the recommended water level

AC Power Consumption - 100 - 240 V AC, 50/60 Hz, Max 2 A  
(Target Tidal Volume Function)

BMC Y series BPAP Systems are Bi-level PAP (Bi-level Positive Airway Pressure) devices intended to provide non-invasive ventilation for patients with Obstructive Sleep Apnea (OSA) and Respiratory Insufficiency. They are intended for adult patients by prescription in the home or hospital/institutional environment. With its Target Tidal Volume function and other excellent comfort features and effective performance, it offers each patient personalized ventilation support.

[Technical Spec-BMC](#)

