

A photograph of a young child with dark hair, wearing a blue long-sleeved shirt and a white face mask with red ties. The child is holding the mask up to their face with both hands, giving a thumbs-up gesture. The background is a blurred outdoor setting.

Buy2Give

Double Qty to charity

BUY2GIVE FACEPROTECTOR 95

€4.95

Buy2Give - Matching Donation.



PEACE 2000 FOUNDATION

Donate any amount for purchasing protectors. Double of the ordered quantity will shipped to charity.

Made of high-quality cotton and non-woven meltblown filter cloth. Multi layer filtration protection that is comfortable to wear. Fits tightly around nose and mouth for maximum protection.

CE and FDA certified FFP2 KN95 Disposable

SKU: DCC95-Buy2Give

Category: [Earthchild](#)

PRODUCT DESCRIPTION

Click for Detailed Description

FFP2 Disposable KN95 FaceProtect

Material: 4ply(non-woven+Hot air cotton+meltblown cloth+non-woven)

Color: White

Standard: China GB2626-2006 KN95 = N95 = FFP2

Application:Personal protection,PM2.5 protection

Filtering effect:95%

Use for hospital,laboratory,food industry,dust free workshop,electronic manufacturers.

Product Features

- 1.Multi layer filter structure protection
- 2.Ear belt / head-mounted
- 3.Comfortable and breathable
- 4.Good protective effect
- 5.Comfortable soft elastic band
- 6.Soft sponge nose pad
- 7.Aluminum adjustable nose clip

General information - Contact us for more details if needed:

Supplied from multiple factories

All masks quality checked and CE certificates verified

CE certificate / Test Certificates / Photos can vary by manufacturer.

Packaging: Variable. Op-Bags 1/4 pcs / Larger bags 10/20/50 pcs.

Typical QTY in carton: 1000 pcs.

Special packing configuration available upon request.

[Bioide-Type95_TechFile220520](#)



