



**GLOVES - POWDER FREE VINYL  
PATIENT EXAMINATION GLOVES  
(100PCS)**

**\$116.82**

100 pcs pack Powder free vinyl patient  
examination gloves

---

**SKU:** CZT-gloves

**Categories:** [Bioide](#), [Gloves](#), [PersonProtect](#)

**PRODUCT DESCRIPTION**

Powder free vinyl patient examination gloves

Product details

type length palm width Thickness of the palm S 240±10mm 85±5mm ≥0.08mm M 95±5mm L  
105±5mm XL 115±5mm

CE and FDA certificates.

Packing: 4000 pcs / box - 10 boxes in carton

Supply Ability: 8 million gloves a month

